



ENDOCLINIC

Advanced Medical Image & Record Management System



Gastrointestinal
Endoscopy



Pulmonary
Endoscopy



Urology
Endoscopy



Laryngoscopy



Laparoscopy



Hysteroscopy
in IVF Centers



Veterinary
Endoscopy

About ENDOClinic

Detailed post-procedure reports are essential for the proper treatment & care of patients.


Managing patient records manually is a cumbersome job with time-consuming reporting, limited ability to display images, manual data backup, etc. In today's digital world, having a software to do all this effortlessly is a given.

Mitra's advanced Medical Image and Records Management System "EndoClinic" offers all the above and more, making your reporting and patient records storage an easy task.


EndoClinic is a versatile and highly efficient management system that helps users achieve better clinical outcomes through accurate documentation & reporting. The software can capture full HD video & unmatched image quality with large-sized images for better visualisation.

Make information work for you.


Evolution of ENDOClinic



Top Features



User-friendly



Full HD Image Viewing






Image Capture




View images in large size



Movie recording, editing and merging




Video Capture




Easy patient record search




Secure database




Import/Export Image & Video



Backup patient data in CD/DVD/Pendrive



Create & edit new templates




Patient data summary in excel




ENDOCLINIC LAN

Collaborate seamlessly with LAN connectivity




ENDOCLINIC Connectivity



ENDOCLINIC Reports


ENDOCLINIC trusted by doctors for well-formatted customised reports

Sample Reports









GASTRO CLINIC

B-226, Okhla Area, New Delhi, India
Ph.: 18001020307



Patient Name: Mrs. Test	Age/Sex: 56 Yrs/F	Ref By: Self
Registration No.: 2	Date: 27-04-2023 14:58:00	Procedure: Colonoscopy
Diagnosis: Test		

COLONOSCOPY REPORT


Findings

Colonoscope passed till terminal ileum
Rectum shows normal mucosa
Sigmoid is normal
Descending Colon is normal
Transverse Colon is normal
Ascending Colon is normal
Cecum is normal
Terminal ileum is normal

Conclusion
Normal Colonoscopy


Dr.
DM Gastroenterologist

Advanced Information Management By EndoClinic(www.endoclinic.in)



PULMONARY CLINIC

B-226, Okhla Area, New Delhi, India
Ph.: 18001020307



Patient Name: Miss. ACD	Age/Sex: 12 Yrs/F	Ref By: Self
Registration No.: 3	Date: 27-04-2023 18:55:00	Procedure: Bronchoscopy
Diagnosis:		


BRONCHOSCOPY REPORT


Findings


Vocal Cords : Normal
Trachea : Normal
Carina : Normal
Right Bronchus : Normal
Left Bronchus : Normal


Conclusion
Normal Bronchoscopy

Dr.
Pulmonologist











Advanced Information Management By EndoClinic(www.endoclinic.in)



UROLOGY CLINIC

B-226, Okhla Area, New Delhi, India
Ph.: 18001020307



Patient Name: Mr. Test A	Age/Sex: 53 Yrs/M	Ref By: Self
Registration No.: 3	Date: 27-04-2023 16:00:00	Procedure: Ureteroscopy
Diagnosis: Test		


CYSTO-URETHROSCOPY REPORT


Findings


Meatus normal
Urethra normal
Bladder neck normal
Trigone normal
RUO normal
LUO normal
Dome normal


Conclusions
Normal Cysto-Urethroscopy



Dr.
Urologist









Advanced Information Management By EndoClinic(www.endoclinic.in)



GI CLINIC

B-226, Okhla Area, New Delhi, India
Ph.: 18001020307



Patient Name: Mrs. S N P	Age/Sex: 45 Yrs/F	Ref By: Self
Registration No.: 7	Date: 27-04-2023 21:57:00	Procedure: Duodenoscopy
Diagnosis:		

ERCP REPORT

Papilla : Normal
Pancreatic Duct : Not Cannulated
Common Bile Duct : 14mm Multiple filling defects seen in the proximal

Technique : CBD cannulated _ Sphincterotomy extended over the wire. Multiple stones around 8-9mm in size removed by extraction balloon. CBD swept multiple times with balloon. Sludge removed. 7Fr x 7cm stent placed.

Conclusions : CBD stent placed

Dr.
Gastroenterologist









Advanced Information Management By EndoClinic(www.endoclinic.in)

Increase productivity with **ENDO**Clinic

Features	EndoClinic 8iHD	EndoClinic 7	EndoClinic 8iHD LAN	EndocCinic 7 LAN
LAN Connectivity	✗	✗	✓	✓
Full HD image & video quality	✓	✗	✓	✗
Multi-user environment	✗	✗	✓	✓
Export/Import images & videos	✓	✓	✓	✓
Well-formatted multi-pages report	✓	✓	✓	✓
Patient data summary in excel	✓	✓	✓	✓
Easy patient record search	✓	✓	✓	✓
Secure database	✓	✓	✓	✓
Intelligent image cropping	✓	✗	✓	✗
Image & Video capture resolution	1920/1080p	720/576p	1920/1080p	720/576p
Compatibility	Windows 10 onwards	Windows 10 onwards	Windows 10 onwards	Windows 10 onwards
Required configuration	i5/i7 Processor, 8 GB RAM, 1 TB Hard disk, USB/PCI slot	i3/i5/i7Processor, 8 GB RAM, 500 GB Hard disk, USB/PCI slot	i5/i7 Processor, 8 GB RAM, 1 TB Hard disk, USB/PCI slot	i3/i5/i7Processor, 8 GB RAM, 500 GB Hard disk, USB/PCI slot

High Definition

ENDOClinic 8iHD

ENDOClinic 8iHD LAN

Standard Definition

ENDOClinic 7

ENDOClinic 7 LAN

Scan QR Code

www.endoclinic.in



**Specifications, Designs and Models are subject to change without any obligation on the part of the manufacturer*

MITRA MEDICAL SERVICES LLP

B-226, Okhla Industrial Area, Phase-1,
New Delhi - 110 020. INDIA

Toll Free No.: 1800 1020 307

E-Mail: endoclinic.services@ottomed.com

Printed - 29 MAY 2023